



Release Notes

Software Update 24.35.01 – vet fluidlab 1

What's new?

On **January 20th, 2025** anvajo will release a new software update for the vet fluidlab 1. This update is designed to improve reliability and ensure smooth performance in your daily use.

What's included in this update?

- **Improved Automatic Slide Recognition:** We've resolved a bug in the slide recognition system to make it more reliable, improving the overall user experience.
- **Spectrometer Firmware Update:** The firmware for the spectrometer has been updated to enhance device stability and maintain optimal performance for the automatic start of the measurement.

This is a maintenance update with no changes on performance or how you use the device – it's simply about ensuring everything works as smoothly as possible.

Please update your device soon to benefit from these improvements. Contact us if you have any questions or if you need help.